

INSTECH Flash Table-Top Autoclave N-Type

An ISO 9001/ISO 13485/CE Company

INSTECH Flash Table-Top Autoclave N-Type conform to the latest world standards.



Model :ISFN-25



Model :ISFN-25

Product Description

- ✓ INSTECH Table-Top Autoclave ideal for Sterile applications in Hospitals, Medical Colleges, Research Centers & Pharma Units.
- ✓ Our manufacturing process entails best-in-class engineering practices like CNC Laser Cutting & bending followed by TIG Welding & Mechanised Finishing.
- ✓ Our products undergo a comprehensive QA regime and come with all safety devices required for a trusted and safe operation.
- ✓ We conduct thorough testing of parts for durability.

Salient Features

- ✓ N-Process.
- ✓ Audio Visual Alarm to Notify Temp./Pressure Deviation, Pressure Sensor ,Power & Door Improper Closed failures.
- ✓ Built-in Open type Water Tank.
- ✓ Automatic Control With LCD Display.
- ✓ Relief Valve in Case of High Pressure.
- ✓ Performs All Required Tests.

Bespoke Customized Capabilities

- FLASH TABLE WITH UNDER SHELF
- SS Tray Shelves 3 Nos.
- Standard USB Ports/ HMI
- Working Pressure (-0.95 - 2.5 bar)
- Sterilization Temp. (121°C- 134°C)
- These all are available Optionallly

MODELS

Model	ISFN-25
Cap.	25 Ltr.
Dinensions (mm)	270Dia X450 Depth
Power Consumption	1950W

INSTECH SYSTEMS

ADD.: Hall No. 5, Plot No. 173/3,
Rajendra Nagar Indl. Area, Ghaziabad, PIN-201007.
Contact : +91 9873069138 +91 8896456000, +91 9812014170
Website: www.cssdtechnologies.com Email: cssd.tecnik.de@gmail.com