

INSTECH ETO Sterilizer conforms to the latest world standards.



Product Description

- ✓ INSTECH ETO Sterilizer is ideal for Sterile applications in Hospitals, Medical Colleges, Research Centres & Pharma Units
- ✓ Our manufacturing process entails best-in-class engineering practices like CNC Laser Cutting & bending followed by TIG Welding & Mechanised Finishing
- ✓ Our products undergo a comprehensive QA regime and come with all safety devices required for trusted and safe operations
- ✓ We conduct thorough testing of parts for durability.

Salient Features

- ✓ Negative Pressure working System.
- ✓ Real time monitoring of the sterilization chamber pressure changes.
- ✓ Effective exhaust of residual Eto.
- ✓ Humidification control system : built in humidity sensor.
- ✓ Safety Door Interlock: during operation, the door doesn't open

Bespoke Customized Capabilities

Automatic Cartridge Puncture Device

Capacity: 50,90,120,160,
200,250,425 Liters.

NO External Input Pressure/Vaccum

Optional: Sealing Station, Baskets

Thermal Printer

56°C HIGH Temp. & 37°C Low Temp.

Models

Model	Dimension	Capacity (Itr)
IS EO 50	300X300X600	50
IS EO 90	300X300x1050	90
IS EO 120	300X300x1350	120
IS EO 160	400X400x1050	160
IS EO 200	400X400x1250	200
IS EO 250	450X450x1250	250
IS EO 425	600X600x1200	425