

INSTECH Flash Table-Top Autoclave B-Type

INSTECH Flash Table-Top Autoclave B-Type conform to the latest world standards.



Model: ISFB-25



Model: ISFB-25

Product Description

- ✓ INSTECH Table-Top Autoclave ideal for Sterile applications in Hospitals, Medical Colleges, Research Centers & Pharma Units.
- ✓ Our manufacturing process entails best-in-class engineering practices like CNC Laser Cutting & bending followed by TIG & Mechanised Finishing
- ✓ Our products undergo a comprehensive QA regime and come with all safety devices required for a trusted and safe operation
- ✓ We conduct thorough testing of parts for durability

Salient Features

- ✓ B-Process, Vacuum Drying Process,
- ✓ Audio Visual Alarm Notifies Temp., Pressure Deviation Pressure Sensor ,Power & Door Improper Closed Failures.
- ✓ Built-in Open type Water Tank.
- ✓ Vacuum drying System Diaphragm VACUUM PUMP.
- ✓ Automatic Control With LCD Display.
- ✓ Relief Valve in Case of High Pressure.
- ✓ Performs B&D test, VLT & Helix test

Bespoke Customized Capabilities

High Efficiency Ultimate Vacuum

Sterilization Temp. (121°C- 134°C)

Standard USB Ports/HMI

Working Pressure(-0.95 - 2.5 bar)

SS Tray Shelves 3 Nos.

These all are available optionally

MODELS

Model	ISFB-25
Cap.	25 Ltr.
Dimensions (mm)	270Dia X 450Depth
Power Consumption	1950W

INSTECH SYSTEMS

ADD.: Hall No. 5, Plot No. 173/3,
Rajendra Nagar Indl. Area, Ghaziabad, PIN-201007.

Contact : +91 9873069138, +91 8896456000, +91 9812014170

Website: www.cssdtechnologies.com Email: cssd.tecnik.de@gmail.com