

**INSTECH Ultrasonic Cleaner conform to the latest world standards.**



**Capacity 30 Ltr.**



**Capacity 20 Ltr.**

## Product Description

- ✓ INSTECH ULTRASONIC CLEANER is ideal for Cleaning applications in Hospitals, Medical Colleges, Research Centers & Pharma Units
- ✓ Use of Best in class manufacturing processes & engineering practices like CNC Laser Cutting & bending, TIG Welding & Mechanised Finishing
- ✓ Our products undergo a comprehensive QA regime and come with all safety devices required for a trusted and safe operation
- ✓ We conduct thorough testing of parts, thus ensuring reliability and durability

## Salient Features

- ✓ Adjustable Digital Timer and Temperature Indicator
- ✓ Highly efficient cleaning can be done using industrial alcohol at High frequency
- ✓ A well built & reliable controller circuit
- ✓ Can be used in the Biochemistry, Microbiology
- ✓ Cleaning to be done using RO water only

## Models

| Models Name | Capacity |
|-------------|----------|
| IS USC-20   | 20 Ltr.  |
| IS USC-30   | 30 Ltr.  |
| IS USC-40   | 40 Ltr.  |
| IS USC-50   | 50 Ltr.  |
| IS USC-80   | 80 Ltr.  |

## Technical Features

|                  |              |
|------------------|--------------|
| Ultrasonic Freq. | 40 KHz       |
| Input Supply     | AC 220-240V  |
| Freq. in Hz      | 50 Hz.       |
| Time Range       | 1-15 min     |
| Ultrasonic Power | 500W         |
| Temp Range       | 30°C to 65°C |
| Heating Power    | 500/1000W    |